EN: Read the precautions for use
FR : Lire les précautions d'emploi
DE: Lesen Sie die Vorsichtsmaßnahmen für die Verwendung
ES: Lea las precauciones de uso.
NL: Lees de voorzorgsmaatregelen voor gebruik
IT: Leggi le precauzioni per l'uso
PT: Leia as precauções de uso
CZ: Přečtěte si opatření k použití
PL: Przeczytaj środki ostrożności dotyczące użytkowania
RO: Citiți precauțiile de utilizare
SV: Läs försiktighetsåtgärderna för användning

PACK-LIPOFILLING

REF: PACK000

FR

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EN

Intended use:

The devices contained in the lipofilling pack are used to perform autologous free fat transfer. The procedure usually involves harvesting adipose tissue with liposuction, preparing the aspirated tissue with various techniques, and injecting the prepared tissue into the receptor sites.

Target Population:

For human use – Do not use on children or pregnant women

Intended Users:

Healthcare professionals only

Contraindications and side effects:

The contraindications and side effects of fat harvesting cannulas are described in the NUP-LT LEAFLET 'collection and injection cannulas for lipotransfer' attached to the kit.

The contraindications and side effects of reinjection cannulas are described in the NUP-MI LEAFLET 'MICRO-CANNULAS AND FILLER INJECTION KITS'.

Instructions for use

The collection and injection of autologous adipose tissue is performed in 1 surgical step under general or local anesthesia.

An anesthetic solution is injected subcutaneously at the intended cannula insertion points.

The incision site is disinfected and a small incision is made using the 14G or 21G needle. Small incisions are made in different places depending on the area to be removed, but are always designed to hide the surgical scar.

The tumescence solution is distributed at the sampling area.

Leave tumescent fluid on for at least 15 minutes.

Perform fat sampling using the CPXXXXXX cannula connected to a syringe (not supplied).

Transfer fat tissue into syringes using the 2.4 mm FLL/FLL connector.

Preparation of the fat alone:

Homogenize the fat graft before injection, connect the 2.4mm connector between the two syringes, and transfer the fat graft from one syringe to the other repeatedly.

The fat graft can be transferred to syringes via the 2.4 mm connector. The micronized fat is ready for injection via the MIXXXX micro cannulas into the area identified by the practitioner after incision with the 21G needle.

Warnings and precautions:

The warnings for harvesting cannulas are described in the package leaflet NUP-LT 'COLLECTION AND INJECTION CANNULAS FOR LIPOTRANSFER' and NUP-MI 'MICRO CANNULAS AND FILLER INJECTION KITS' attached to the kit.

The warnings and precautions of other devices are described below:

PRECAUTIONS OF THE DEVICES CONTAINED IN THE KIT:

Read the instructions and precautions on each device before use.

Precautions for sterile single-use hypodermic needles, reference 8300013707 and 4710008040:

- If the packaging is damaged or opened, the product should not be used due to the potential alteration of sterility conditions.
- Connect the hypodermic needle to a suitable medical device (e.g. a syringe).
- Make sure the needle is securely attached, then carefully remove the protector.
- Do not recap. Dispose of the needle in a sharps container.
- Once used, do not reuse or resterilize it.
-

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