EN: Read the precautions for use
FR : Lire les précautions d'emploi
DE: Lesen Sie die Vorsichtsmaßnahmen für die Verwendung
ES: Lea las precauciones de uso.
NL: Lees de voorzorgsmaatregelen voor gebruik
IT: Leggi le precauzioni per l'uso
PT: Leia as precauções de uso
CZ: Přečtěte si opatření k použití
PL: Przeczytaj środki ostrożności dotyczące użytkowania
RO: Citiți precauțiile de utilizare
SV: Läs försiktighetsåtgärderna för användning

PACK-LIPOFILLING

REF: PACK-MELT-1

REF: PACK-MALT-1

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[FABRICANT]

Shanghai Kindly Medical Instruments Co, Ltd

Adresse : No. 925 Jinyuan y : No. 925 Jinyuan yi Road, 201803 Shanghai, P.R.C. Tel : +86-(0)21-59140056

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Shanghai International Holding Corp, GmbH (Europe)

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Tel : 0049-40-2513175 MADE IN CHINA

Intended Use:

The devices contained in the lipofilling pack are used to perform autologous free fat transfer. The procedure generally involves harvesting adipose tissue via liposuction, preparing the aspirated tissue using various techniques, and injecting the prepared tissue into the recipient sites.

Target Population:

For human use only – Do not use on children or pregnant women.

Intended Users:

For healthcare professionals only.

Contraindications and Side Effects:

The contraindications and side effects of infiltration and harvesting cannulas are described in the NUP-LT leaflet '*HARVESTING AND INJECTION CANNULAS FOR FAT TRANSFER*' included in the kit. The contraindications and side effects of reinjection cannulas are described in the NUP-MI leaflet '*MICRO-CANNULAS AND FILLER INJECTION KITS*'.

Instructions for Use:

The devices included in the lipofilling pack are intended for performing an autologous free fat transfer. This procedure typically involves three steps: harvesting adipose tissue via liposuction, preparing the aspirated tissue using various techniques, and injecting the prepared tissue into the recipient sites.

1. Harvesting Adipose Tissue

- The collection and injection of autologous adipose tissue are performed in a single surgical procedure under general or local anesthesia.
- The incision site is disinfected, and a small incision is made using a **14G needle (8300013707 / AF1440)**.
- Small incisions may be made at various points depending on the harvesting area while being positioned to minimize visible scarring.
- The tumescent solution is infiltrated into the harvesting area using the **CI21150 infiltration cannula**.
 - Infiltration can be performed directly with the **30 ml syringe (ZZSQ30)** provided in the kit.
 - Alternatively, the cannula can be connected to a bag containing the infiltration solution via a stopcock; in this case, the syringe serves as a transfer device. **(The infiltration solution and stopcock are not included in the kit.)**
- Allow the tumescent solution to act for **at least 15 minutes** before proceeding with the harvesting.
- Fat harvesting is performed using the **CP141506HS cannula**, connected to a **30 ml syringe (ZZSQ30)**.

2. Preparation of Adipose Tissue

- The harvested adipose tissue is transferred into **10 ml syringes (BD® Luer-Lock™ Tip 300912 or 305959) or 1 ml syringes (309628)** using an **FLL/FLL connector**:
 - **FD0000000-LLF12 (1.2 mm)**
 - **FD0000000-LLF14 (1.4 mm)**
 - **FD0000000-LLF24 (2.4 mm)**
- **Sedimentation or Centrifugation:**
 - Allow the graft to sediment using **LL caps FD0001711-MC**.
 - Alternatively, centrifuge the sample after transferring it into a syringe or a tube dedicated to centrifugation **(not included in the kit)**.

- Once centrifugation is complete, remove the syringe while keeping it in a vertical position to avoid mixing the layers.

- **Phase Separation:**

- After centrifugation, three layers can be observed:
 1. The top layer contains lipids.
 2. The middle layer consists of adipose tissue.
 3. The bottom layer contains blood, tissue fluid, and the infiltration solution.
- Only the middle layer of adipose tissue should be retained.
- Remove the top layer using an absorbent material or by tilting the syringe.
- Eliminate the bottom layer by carefully draining the syringe from the lower end.

3. Final Preparation and Injection

- **Fat Graft Homogenization:**

- Before injection, connect a **2.4 mm connector (FD0000000-LLF24)** between two **10 ml syringes (BD® Luer-Lock™ Tip 300912 or 305959) or 1 ml syringes (309628)** and transfer the graft back and forth multiple times to achieve optimal homogenization.

- **Micronized Fat Preparation:**

- The fat graft can be transferred into **1 ml syringes (309628)** using a **2.4 mm connector (FD0000000-LLF24) or 1.4 mm connector (FD0000000-LLF14)**.
- If further particle size reduction is needed, use a **1.4 mm connector (FD0000000-LLF14) or 1.2 mm connector (FD0000000-LLF12)**.

- **Injection:**

- The micronized fat is ready for injection using the **CF201590 or CF14790 injection cannula** into the target area as determined by the practitioner, after an incision made with a **14G needle (8300013707 / AF1440)**.

Warnings and Precautions:

The warnings related to infiltration, harvesting, and reinjection cannulas are described in the NUP-LT leaflet 'HARVESTING AND INJECTION CANNULAS FOR FAT TRANSFER' and the NUP-MI leaflet 'MICRO-CANNULAS AND FILLER INJECTION KITS' included in the kit. The warnings and precautions for other devices are detailed below:

PRECAUTIONS FOR THE DEVICES CONTAINED IN THE KIT:

Read the instructions and precautions for each device before use.

Precautions for the BD 10 ml Luer-Lok™ Tip Syringe, reference 305959 or 300912:

Warning : To prevent the transmission of HIV (AIDS), HBV (hepatitis), and other infectious diseases from accidental needle sticks, used needles should not be recapped or removed unless absolutely necessary or required by a specific medical procedure. OSHA (USA) standards require that needle recapping be performed using a one-handed technique.

Dispose of all needles and syringes in accordance with the facility's regulations.

Pyrogen-free. Do not use if the individual packaging is damaged. Do not use for oral injections.

Caution: Reuse may cause infection or other adverse effects.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

This product is not made with natural rubber latex.

Caution: Pharmaceuticals prepared in these syringes should be administered as soon as possible.

Precautions for sterile single-use hypodermic needles, reference 8300013707 :

If the packaging is damaged or open, the product must not be used due to potential sterility compromise.

Connect the hypodermic needle to an appropriate medical device (e.g., a syringe).

Ensure the needle is securely attached, then carefully remove the protector.

Do not recap the needle. Dispose of it in a sharps container.

Do not reuse or resterilize after use.

Instructions for Use of the Shanghai Kindly Enterprise Development Group 30ml Syringe, Reference ZZSQ30

[PRODUCT NAME] Needle-Free Syringes

Needle-free syringes consist of a mandrel, a mandrel cap, a sleeve, a handle, a piston, and a connector. The raw materials used are PC, rubber, PE, and ABS.

[PRODUCT TYPE] 1ml, 3ml, 6ml, 10ml, 20ml, 30ml

[APPLICATION FIELD]

Sterile, single-use needle-free syringes, sterilized by EO, are mainly used to assist blood vessels in manifesting by pushing the liquid into the vascular channel during surgical procedures.

[INSTRUCTIONS FOR USE]

1. Open the packaging and remove the product.
2. Inject the liquid and expel the air.
3. Connect the device to other devices, then push the liquid into the vascular channel to help manifest the blood vessels.

[WARNINGS]

1. Reusing single-use devices may pose a potential risk to the patient or user. It may lead to contamination and/or a loss of functionality. Contamination and/or loss of functionality of the device may cause injury, illness, or death to the patient.
2. Medical waste must be disposed of in accordance with local laws and regulations.

[PRECAUTIONS]

1. If the individual packaging is damaged, do not use it!
2. The product is for single use only. Dispose of it after use; do not reuse!
3. If the product is expired, do not use it!
4. Do not alter the structure of the product during use!

[CONTRAINDICATIONS]

There are no absolute contraindications.

[STORAGE]

The product should be stored away from sunlight and rain, in a cool and dry place.

[SYMBOLS AND LABELING]

Do not reuse

Sterilized by EO

Caution, refer to instructions for use.

[SHELF LIFE] Three years.

[MANUFACTURER]

Shanghai Kindly Medical Instruments Co., Ltd

Address: No. 925 Jinyuan Yi Road, 201803 Shanghai, P.R.C.

Tel: +86-(0)21-59140056

[EUROPEAN REPRESENTATIVE]

Shanghai International Holding Corp, GmbH (Europe)

Address: Eiffestraße 80, 20537 Hamburg, Germany

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