EN: Read the precautions for use
FR : Lire les précautions d'emploi
DE: Lesen Sie die Vorsichtsmaßnahmen für die Verwendung
ES: Lea las precauciones de uso.
NL: Lees de voorzorgsmaatregelen voor gebruik
IT: Leggi le precauzioni per l'uso
PT: Leia as precauções de uso
CZ: Přečtěte si opatření k použití
PL: Przeczytaj środki ostrożności dotyczące użytkowania
RO: Citiți precauțiile de utilizare
SV: Läs försiktighetsåtgärderna för användning

PACK-LIPOFILLING

REF: PACK-MILT

REF: PACK-NALT

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Intended Use:

The devices contained in the lipofilling pack are used to perform autologous free fat transfer. The procedure generally involves harvesting adipose tissue via liposuction, preparing the aspirated tissue using various techniques, and injecting the prepared tissue into the recipient sites.

Target Population:

For human use only – Do not use on children or pregnant women.

Intended Users:

For healthcare professionals only.

Contraindications and Side Effects:

The contraindications and side effects of infiltration and harvesting cannulas are described in the NUP-LT leaflet '*HARVESTING AND INJECTION CANNULAS FOR FAT TRANSFER*' included in the kit. The contraindications and side effects of reinjection cannulas are described in the NUP-MI leaflet '*MICRO-CANNULAS AND FILLER INJECTION KITS*'.

Instructions for Use:

The collection and injection of autologous adipose tissue are performed in a single surgical step under general or local anesthesia.

An anesthetic solution is injected subcutaneously at the planned cannula insertion points.

The incision site is disinfected, and a small incision is made using a 14G needle. Small incisions are made at various points depending on the harvesting area, but they are always designed to minimize visible scarring. The tumescent solution is distributed in the harvesting area using the infiltration cannula CXXXXX. Allow the tumescent fluid to act for at least 15 minutes.

Harvest the fat using the CPXXXXX cannula connected to a 10 ml syringe. Transfer the adipose tissue into 1 ml syringes using the FLL/FLL connector (2.4 mm, 1.4 mm, or 1.2 mm). Wait for sedimentation using the FD0001711-MC caps or centrifuge the fat sample after transferring it into a syringe or tube intended for centrifugation (not included in the kit). Remove the syringe or tube carefully, keeping it in an upright position to avoid mixing.

After centrifugation, three layers can be observed:

- The first layer contains lipids,
- The second layer consists of adipose tissue,
- The third layer contains blood, tissue fluid, and infiltration solution.

Only retain the adipose tissue by removing the first layer with an absorbent material or by tilting the syringe and eliminating the third layer from the bottom.

Preparation of purified fat:

To homogenize the fat graft before injection, connect the 2.4 mm, 1.4 mm, or 1.2 mm connector between two 10 ml or 1 ml syringes and transfer the fat graft back and forth multiple times. The fat graft can be transferred into 1 ml syringes using the 2.4 mm or 1.4 mm connector. If further particle size reduction is needed, the 1.4 mm or 1.2 mm connector can be used. The micronized fat is ready for injection via micro-cannulas MIXXXX in the area identified by the practitioner after incision with a 21G needle.

Warnings and Precautions:

The warnings related to infiltration, harvesting, and reinjection cannulas are described in the NUP-LT leaflet '*HARVESTING AND INJECTION CANNULAS FOR FAT TRANSFER*' and the NUP-MI leaflet '*MICRO-CANNULAS AND FILLER INJECTION KITS*' included in the kit.

The warnings and precautions for other devices are detailed below:

PRECAUTIONS FOR THE DEVICES CONTAINED IN THE KIT:

Read the instructions and precautions for each device before use.

Precautions for the BD 10 ml Luer-Lok™ Tip Syringe, reference 305959 or 300912:

Warning : To prevent the transmission of HIV (AIDS), HBV (hepatitis), and other infectious diseases from accidental needle sticks, used needles should not be recapped or removed unless absolutely necessary or required by a specific medical procedure. OSHA (USA) standards require that needle recapping be performed using a one-handed technique.

Dispose of all needles and syringes in accordance with the facility's regulations.

Pyrogen-free. Do not use if the individual packaging is damaged. Do not use for oral injections.

Caution: Reuse may cause infection or other adverse effects.

Caution: Federal law (USA) restricts this device to sale by or on the order of a physician.

This product is not made with natural rubber latex.

Caution: Pharmaceuticals prepared in these syringes should be administered as soon as possible.

Precautions for sterile single-use hypodermic needles, reference 8300013707 :

If the packaging is damaged or open, the product must not be used due to potential sterility compromise.

Connect the hypodermic needle to an appropriate medical device (e.g., a syringe).

Ensure the needle is securely attached, then carefully remove the protector.

Do not recap the needle. Dispose of it in a sharps container.

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